



# Real-Time Experimental Control

(REC™) for live cell fluorescence

reporter applications.

Application Note

**AUTOMATED MONITORING OF THE REPORTER ACTIVITY OF HFOB1.19 CELLS  
FOLLOWING INDUCTION OF OSTEOGENIC DIFFERENTIATION BY TEMPERATURE SHIFT**





Parameter	Setting
Temperature control	On
Target temperature	34 °C
CO <sub>2</sub> control	On
Target conc.	5 %
Plate	TEC96ft_cell
Mode	Kinetic
Kinetic cycles	41
Interval time (hh:mm:ss)	04:00:00
Measurement mode	Fluorescence imaging
Application	User defined
Channels	Bright field, Green
Objective	4x
Pattern	User defined
<b>Bright field channel</b>	
Focus offset	-15 µm
Sensitivity	50 %
Object size	8-30 µm width, 8-30 µm length
Fill holes, with size <x µm <sup>2</sup>	0 µm
<b>Green channel</b>	
LED intensity	100 %
Focus offset	0 µm
Exposure time	300 ms
Sensitivity	50 %
Object size	8-30 µm width, 8-30 µm length
Analysis type	Segmentation
<b>Real-time control</b>	
Mode	Condition
Type	Start at value
Executed once	
Input data	Label 1 Confluence
Reference	A2
Operation	>80 %
Mode	Condition
Type	Start at cycle
Operation	At cycle, 10, 20, 30, 40, 50
User intervention	Check humidity and medium change

Table 1: Software protocol to assess reporter activity on the induction of osteogenic differentiation by temperature shift.

## RESULTS

hFOB1.19 cells cultured inside the Spark Cyto reader showed a steady and homogeneous growth over the plate (CV = ≤10 %) from 20 to 80 % confluence. Cell growth continued, showing identical dynamics, following an automated temperature shift from 34 to 39 °C at 80 % confluence, finally reaching 99 % confluence after 183 hours (see Figure 1). The temperature shift led to the expression of eGFP in the cells, increasing the green fluorescent signal significantly in all samples (see Figure 2).




Figure 1: Growth of hFOB1.19 cells with a kinetic condition at 80 % confluence triggering a temperature increase to 39 °C.




Figure 2: Induction of eGFP expression in hFOB cells by temperature shift from 34 °C (cycle 30) to 39 °C (cycle 41).

## CONCLUSIONS

The results acquired in this study demonstrate the capability of the Spark Cyto to conduct long-term live cell assays, such as live cell fluorescence reporter applications, for up to seven days. This is due to the full set of integrated environmental control features, including temperature regulation, gas control and evaporation protection. Special software features, such as kinetic conditioning, help to minimize hands-on time by automating critical steps of the workflow, for example the induction of temperature shifts at a certain confluence level.



## ABBREVIATIONS

DMEM	Dulbecco's modified Eagle's medium
eGFP	enhanced green fluorescent protein
FBS	fetal bovine serum
hFOB1.19	human fetal osteoblast 1.19
LED	light-emitting diodes
MEM	modified Eagle's medium
REC	Real-Time Experimental Control

### About the author

*Dr. Christian Oberdanner is Senior Application Specialist at Tecan Austria. He studied molecular biology at the University of Salzburg with a strong focus on cell- and tumor biology. Christian started to work for Tecan Austria as external scientific consultant in 2005 and permanently joined the company in 2006. Since then he held several roles as Application Scientist, Application Specialist and Product Manager in research and development as well as in the sales and marketing department. Christian's priority within Tecan are multimode microplate reader applications and cell imaging.*

### For Research Use Only. Not for use in diagnostic procedures.

.....  
**Australia** +61 3 9647 4100 **Austria** +43 62 46 89 330 **Belgium** +32 15 42 13 19 **China** +86 21 220 63 206 **France** +33 4 72 76 04 80 **Germany** +49 79 51 94 170  
**Italy** +39 02 92 44 790 **Japan** +81 44 556 73 11 **Netherlands** +31 18 34 48 17 4 **Nordic** +46 8 750 39 40 **Singapore** +65 644 41 886 **Spain** +34 93 595 25 31  
**Switzerland** +41 44 922 89 22 **UK** +44 118 9300 300 **USA** +1 919 361 5200 **Other countries** +41 44 922 81 11  
 .....

Tecan has not independently validated this method with all possible sample types or analytical uses and is providing this example method as a convenience to users. The user must independently evaluate and validate: (a) the suitability of the method for their use, (b) their ability to process samples of their choosing following the method; and (c) their ability to proficiently perform the method in their facility with their personnel.

Tecan Group Ltd. makes every effort to include accurate and up-to-date information within this publication; however, it is possible that omissions or errors might have occurred. Tecan Group Ltd. cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice. All mentioned trademarks are protected by law. For technical details and detailed procedures of the specifications provided in this document please contact your Tecan representative. This brochure may contain reference to applications and products which are not available in all markets. Please check with your local sales representative.

Tecan, Spark Cyto and SparkControl are registered trademarks or trademarks of Tecan Group Ltd., Männedorf, Switzerland.

© 2019 Tecan Trading AG, Switzerland, all rights reserved. For disclaimer and trademarks please visit [www.tecan.com](http://www.tecan.com)

[www.tecan.com](http://www.tecan.com)

