



TEMPO® parameters

TEMPO parameters	NEXT DAY SOLUTION	AOAC	ISO
41113 AC Aerobic mesophilic total flora	✓	✓	✓
80002 STA Staphylococcus (coag+)	✓	✓	✓
80003 EB Enterobacteriaceae	✓	✓	✓
80004 EC Escherichia coli	✓	✓	✓
80006 TC Coliforms	✓		✓
80044 CC Coliforms	✓	✓	
80406 BC Bacillus cereus	✓	✓	✓
80001 YM Yeast & Mold	72-76 hr	✓	
80071 LAB Lactic acid bacteria	40-48 hr		
421509 CAM Campylobacter	44-48 hr		Coming soon!

TEMPO®

FOOD QUALITY AUTOMATION

Because you count




GET THE RIGHT SOLUTION

SAMPLE PREPARATION



PRODUCTIVITY



SAFETY



IDENTIFICATION



©2018 bioMérieux, Inc. • BIOMÉRIEUX and the BIOMÉRIEUX logo, TEMPO, GENE UP, VIDAS, SMASHER, DILUMAT, CHEMUNEX DCOUNT and VITEK 2 are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies • Any other name or trademark is the property of its respective owner • bioMérieux S.A. RCS Lyon 673 620 399 • theria • RCS Lyon B 398 160 242 • Printed in USA • Patents: www.biomerieux-usa.com/patents • PRN 18-0044-00


SIMPLIFY YOUR LIFE

ONLY 1 CARD


ONLY 3 STEPS

as simple as 1, 2, 3!

- 1 Automated sample preparation



- 2 Incubation in your standard incubators



- 3 Automated reading and recording



HIGHLY FLEXIBLE!


- MANAGEMENT OF SURGE TESTING NEEDS
- FAST START FOR "LAST MINUTE" SAMPLES
- QUICK TECHNICAL TRAINING

"WITH TEMPO, I CAN TEST ON DEMAND"

Sophie CORNU,
Häagen Dazs, General Mills

REDUCE

- WASTE
Protect the environment and be lean
- INCUBATOR SPACE
- SUBJECTIVITY



RESULTS YOU CAN COUNT ON

GLOBALLY RECOGNIZED

- ISO & AOAC VALIDATED
- NO CROSS CONTAMINATION: SECURED PROTOCOL
- LONG SHELF LIFE REAGENTS

DATA INTEGRITY

KEEP YOUR DATA WITHIN COMPLIANCE STANDARDS

- FULL TRACEABILITY - INTEGRATED CONTROL PLAN & BARCODE
- SECURED DATA - FROM SAMPLE ID TO AUTOMATED RESULTS
- NO TRANSCRIPTION ERRORS

— bioMérieux —
CONNECT-UP™

MINIMIZE ADMINISTRATIVE TASKS,
GAIN VALUE-ADDED ACTIONS!